

Water peeling II Clean pores Nutrition introduction Manual▶	INSTRUCTION▶ Cleansing mode 1.Long press 2s power button to power on 2.Press once to select the Cleansing mode Cleansing mode with color light function. Hold the conductive strips on both sides to cleansing face with positive ion. It is recommended to use cleanser first before cleansing the face. 3.Apply the skin care products to the face 4.As shown in Figure 2, to face lifting Figure 1 Figure 2 Moisturizing mode 1.Long press 2s power button to power on 2.Press twice to select Moisturizing mode Moisturizing import mode with color light function Hold the conductive strips on both sides and enter the negative ion mode It is recommended to use a mask first, before using Moisturizing mode to import nutrients. 3. Apply the skin care products to the face 4. As shown in Figure 2, to face lifting Figure 1 Figure 2
---	--

-1-