

Water peeling II

Clean pores Nutrition introduction

Manual ▶

INSTRUCTION ▶

Cleansing mode

1. Long press 2s power button to power on
2. Press once to select the Cleansing mode
3. Apply the skin care products to the face
4. As shown in Figure 2, to face lifting

Moisturizing mode

1. Long press 2s power button to power on
2. Press twice to select Moisturizing mode
3. Apply the skin care products to the face
4. As shown in Figure 2, to face lifting

INSTRUCTION ▶

repair mode


1. Long press 2s power button to power on
2. Press three times to select repair mode
3. Apply the skin care products to the face
4. As shown in Figure 2, to face lifting

MAINTENANCE AND STORAGE ▶

1. Make sure that the unit's power is off.
2. Use a wet tissue or alcohol on the machine head to carefully wiping. Do not wipe this product with corrosive chemical industrial products such as thinner solvents.
3. After clean the machine head, place it in a dry area to safekeeping. Do not store in a damp place such as a bathroom. Do not store in a place where children can touch it.

CHARGING ▶

1. USB charging cable to connect USB interface
Please know the product is arrived without fully charged. Please charge before use.
2. Plug to the cell phone adapter, power bank or computer USB interface.
The LED light flashes continuously during charging, and the LED light is keep turn on when it is fully charged.
Please confirm the status of the lithium battery before charging. Before using for the first time, please remove the sealing film of the skin.



CHARGING PRECAUTIONS ▶

When the battery is low, product will emit a "beep" sound, and then turn off. Note that the remaining amount of the battery is insufficient. Please charge the product immediately. It cost about 2.5-3 hour to be fully charged.

STRUCTURE ▶



FUNCTION ▶

1. Deep cleansing and removing dead skin.
2. Fade pigmentation spots.
3. Anti-wrinkle, enhance skin elasticity.
4. Promote blood circulation and lymph fluid metabolism.
5. Enhance the regeneration ability and permeability of skin cells.
6. face massage , nutrient import, variegated export.

PRODUCT PARAMETERS ▶

Name: Water peeling II
 Input voltage: DC 5V Charging time: about 3H
 Power: ≤3W Instrument net weight: 122g
 Size: 172*54*26mm Working time: 15 minutes
 Battery: 500mA Use time : about a week

PRODUCT FEATURES ▶

Ultrasound is a high-frequency mechanical vibration wave that is widely used in medical and cosmetic fields. The Mechanical waves have mechanical action, warming action and catalysis, which can be used in low-energy high-frequency vibration. The ultrasonic can produce a gentle resonance with skin that enhances cell viability and interstitial tissue Permeability, dredge sweat glands, shrink pores.

Ultrasonic skin scrubber uses Ultrasonic technology which frequency are more than 20,000 times per second. The vibration causes intense friction between the molecules, which can penetrate deep vibration into the skin and make the molecules sharply Friction, deep temperature rise, cell metabolism is strong, micro circulation is accelerated, and metabolite reflux is promoted. Through The conversion of the ion stream, the permeability of the cell membrane is enhanced, and the semi-permeable membrane produces an intracellular slurry microfluid. The fatty acid is catabolized into a fine volume of CO2 and H2O fat. At the same time, it is matched work with some skin care products, beauty effects are more significant.

Ultrasonic skin scrubber through high-frequency vibration to improve microcirculation during fat movement, burning and cleaving the fat cells which converted into metabolites and excreted to achieve deep cleansing. Pick with a special shovel type probe, making ultrasonic vibration effectively enter under skin 4-7 mm to deep cleansing. The liquid becomes finer and atomized, and these dirt and excess sebum can be emulsified to block in the pores if without cleansing. Under ultrasonic vibration, the cosmetic residue, emulsified dirt, and with the splash of water splashing, this cleaning effect is faster and more effective. This is like an ultrasonic cleaner that cleans precision machinery. This phenomenon is also will happen to the skin. Thereby removing the stratum corneum, fat granules, pimples, etc., so that your skin care products are no longer wasted, allowing the skin to absorb better.

WARRANTY CARD ▶

1. Free maintain for the breakdown caused by product quality problems.
2. One year guarantee from the date of purchase.
3. Free to change broken parts, freight is collected by customers.
4. If breakdown caused by the force majeure, free maintain won't be given and freight shall be collected by customers.
5. Please fill the guarantee clearly, without any revision. Please take good care of guarantee. If lost, won't be re-provided.

保証書	
品名	Water peeling II
お名前	
ご住所 〒	
TEL	
ご購入日	年 月 日 保証期間/ご購入日より1年間
販売店印	

発売元 株式会社 エヴァンジェール 製造販売元 株式会社 アグレックス
 ■本社 ■本社
 〒321-0164 〒261-0023
 栃木県宇都宮市双葉1-13-50 千葉県千葉市美浜区中瀬1-3 MTG B棟 9F
 TEL 0120-169-119 TEL 043-350-0111 FAX 043-350-0114